

\*\* February 2024 (12th edition)  
 \* October 2023 (11th edition)

Program 1: Program for Disease Diagnosis  
 Controlled medical device: Program for general-purpose diagnostic imaging device workstation (70030012)

## Medical imaging software Holoeyes MD

### Configuration, structure, principle, etc.


#### Configuration and structure

This program processes imaging information obtained from diagnostic imaging devices and makes it available for healthcare applications. Automatic diagnosis is not performed. In addition to the standard function of image display, it has the optional function of three-dimensional (3D) image processing. The basic functions include an image- and information-processing function (zooming in and out, and movement), and an image-display function.

The program is used after installation on general-purpose IT equipment, etc. It is provided as a download.

#### Operation Principles

The program receives images from a diagnostic imaging device or a server storing such imaging information, and displays them for healthcare purposes.



#### Additional functions

Item	Function
* 3D image-processing function	A function that processes a sequence of imaging data as a 3D image for display. Including virtual endoscopy, etc.
General image measurement function	A basic distance measurement function that uses the pixel values and position information of imaging data.

#### Not applicable to medical devices

Item	Function
Program for educational use	Remote conferencing function that shares 3D images.

### Purpose or effects of use

The program processes human body imaging information from diagnostic imaging devices for healthcare applications. This program does not have an automatic diagnosis function.

### Methods of use, etc.

#### Methods of use

1. Operating environment and preliminary preparation  
 This program is installed on general-purpose IT equipment, etc., that meets the below specifications. Installation is by the method specified by the marketing authorization holder.

① General-purpose IT equipment (recommended environment)

CPU: Intel Core i5-7500 or higher  
 Memory: 8 GB or more  
 GPU: NVidia GTX 1060 or higher  
 HDD: 128 GB or more  
 OS: Windows 10 or higher  
 Safety: Compliant with IEC 60950-1 (JIS C 6950-1) and IEC

62368-1 (JIS C 62368-1)

EMC: Compliant with CISPR32 (or CISPR22) / CISPR24 (or VCCI)

② Microsoft HoloLens

CPU: Intel 32-bit architecture  
 HPU: Custom-built Microsoft Holographic Processing Unit (HPU 1.0)  
 RAM: 2.0 GB  
 Storage: 64 GB Flash  
 Display: 2 HD 16:9 light engine  
 OS: Windows 10 or higher  
 Safety: Compliant with IEC 60950-1 (JIS C 6950-1)  
 EMC: Compliant with CISPR32, FCC/CISPR35 and CISPR24

③ Microsoft HoloLens 2

SoC: Qualcomm Snapdragon 850 computing platform  
 HPU: Second-generation custom-made holographic mounting device  
 RAM: 4.0 GB  
 Storage: 64 GB, UFS 2.1  
 Display: 2k 3:2 light engine  
 OS: Windows Holographic Operating System  
 Safety: Compliant with IEC 62368-1 (JIS C 62368-1)  
 EMC: Compliant with CISPR32, FCC/CISPR35 and CISPR24

④ Magic Leap One

CPU: Parker SOC CPU (NVIDIA<sup>®</sup>) CPU: Denver 2.0  
 RAM: 8.0 GB  
 Storage: 128 GB  
 OS: Lumin OS 0.98  
 Safety: Compliant with IEC 62368-1 (JIS C 62368-1)  
 EMC: Compliant with CISPR32/CISPR24

⑤ Meta Quest

\*\*SoC: Qualcomm<sup>®</sup> Snapdragon 835  
 Graphic: Adreno 540  
 RAM: 4.0 GB  
 Storage: 64/128 GB  
 Display: PenTile OLED 1440 × 1600 per eye @ 72 Hz  
 OS: Android 7 or higher  
 Safety: Compliant with IEC 62368-1 (JIS C 62368-1)  
 EMC: Compliant with CISPR32/CISPR35 and CISPR24

⑥ Meta Quest 2

\*\*SoC: Qualcomm<sup>®</sup> Snapdragon XR2 Platform  
 Graphic: Adreno 650  
 RAM: 6.0 GB  
 Storage: 64/128/256 GB  
 Display: Single high-speed switch LCD 1832 × 1920 per eye @ 72 Hz  
 OS: Android 10 or higher  
 Safety: Compliant with IEC 62368-1 (JIS C 62368-1)  
 EMC: Compliant with CISPR32/CISPR35 and CISPR24

⑦ Meta Quest Pro

\*\*SoC: Qualcomm<sup>®</sup> Snapdragon XR2 Platform  
 Graphic: Adreno 650  
 RAM: 12 GB  
 Storage: 256 GB  
 Display: LCD 1800 × 1920 per eye @ 90 Hz  
 OS: Android 10 or higher  
 Safety: Compliant with IEC 62368-1 (JIS C 62368-1)

Be sure to check the instruction manual.

EMC: Compliant with CISPR32/CISPR35 and CISPR24

**\*\*⑧ Meta Quest 3**

SoC: Qualcomm® Snapdragon XR2 Gen2

RAM: 8 GB

Storage: 128 / 512 GB

Display: Single high speed switch LCD 2064 × 2208 per eye @ 120 Hz

OS: Android 10 or higher

Safety: Compliant with IEC 62368-1 (JIS C 62368-1)

EMC: Compliant with CISPR32/CISPR35

**\* Specifications of external output device**

Minimum 1280 × 720 pixel resolution

32-bit or higher color display

※ When general-purpose IT equipment, etc., with this program downloaded and installed, is used in an environment with patients present, use the internal power supply for the general-purpose IT equipment, etc., in order to ensure safety (JIS T 0601-1, 16.6 Leakage current). Do not use the equipment when connected to an external power supply.

**2. Preparation for use**

- (1) Turn on the power of the device where the program will be installed.
- (2) Start the program.

**3. Operation**

- (1) Obtain imaging data.
- (2) Select a function.
- (3) Display images, etc.

**4. Termination**

- (1) Terminate the program by clicking the termination icon on the screen or selecting the termination option from the menu.
- (2) Turn off the power if necessary.

**Cautions with use**

**Important cautions**

1. Confirm that the 3D model is the appropriate one for the relevant patient before downloading it.
2. Confirm in advance that the images and data are those intended for the relevant patient.
3. Confirm that the display on the external output device screen is not distorted or missing.
4. If your PC freezes for some reason when the 3D model is being displayed, download the 3D model again after verifying that your PC has normal operation.
5. Confirm that the 3D model is not displayed with mirror reversal.
6. Confirm that the 3D model is not displayed smaller or larger than the real thing.
7. If you are too close to a magnified 3D model, the display may not function properly.
8. The educational-use program function is a non-medical function. Do not use it for healthcare.
9. The user is responsible for correcting the images and data.
10. This program should be used by healthcare professionals who are familiar with the collection and processing of imaging information, and display of such images, by each diagnostic imaging device.
11. Obtaining accurate and optimal imaging and data requires proper implementation of all procedures from image acquisition to display, including correct position information, and patient identification. It is the user's responsibility to determine whether the images and data generated are suitable for the intended use.
12. The results displayed are calculated from the images, data, and supplementary information. These may differ

from the actual results due to various factors, such as inappropriate use.

13. The results of measurement by this program are calculated from the images and supplementary information. Therefore, do not make a definitive diagnosis solely on the basis of the measurements obtained with this product, but use them for reference.
14. Do not use this program for primary diagnosis.
15. Do not use this program other than for the purposes of use.
16. Do not use it in a high-concentration oxygen atmosphere, in an explosive atmosphere such as one containing flammable anesthetic gas, or in the presence of inflammable substances.
17. Before using the product, be sure to read the instruction manual for the general-purpose IT equipment, etc.
- \* 18. This program should be downloaded to the general-purpose IT equipment, etc., specified in "Methods of use". When used in combination with other devices, confirm that the device is functioning properly before use in healthcare.
- \*\*19. General-purpose IT equipment on which this program is installed should be stored in a physically controlled location with a key or other means to prevent unauthorized use.
- \*\*20. To prevent unauthorized use, use the account management functions provided with the operating system of general-purpose IT equipment, etc., whenever possible.
- \*\*21. When using this program, please use it in a secure network environment such as in a hospital where security is ensured.
- \*\*22. When connecting to a USB port of a general-purpose IT equipment or other device, make sure that the security of the device to which it is connected is secured before use.

**Contact Information for the Marketing Authorization Holder and Manufacturer**

Holoeyes Inc.

Mourin Building 303, 2-17-3 Minami-Aoyama, Minato-ku, Tokyo

Tel.: 050-1744-9255

\*\* ※ If you encounter a problem that is related to cyber security, please contact us through the contact page on the [Holoeyes MD service website](#).

Be sure to check the instruction manual.